

Plano

Diameter	80mm
Base Curves	6 and 8
Center Thickness	2.0mm


Semi Finished Single Vision

Diameter	70mm
Base Curves	0.50, 2.00, 4.25, 5.25, 6.25, 8.25
Power Range	+7.00 to -10.00 sph, cyl to -4.00

Semi Finished Bifocal



Flat Top 28

Diameter	70mm
Base Curves	2.00, 4.25, 5.25, 6.25, 8.25
Additions	0.75D to 3.50D in 0.25D steps
Power Range	+7.00 to -7.25 sph, cyl to -4.00
Segment Size	28 mm
Segment Placement	5mm in, 5mm down



Progressive

All Purpose F	
Diameter	75mm Decentered 2.5mm
Base Curves	1.50, 3.25, 5.25, 7.25
Additions	1.00D to 3.50D in 0.25D steps
Power Range	+6.50 to -8.00 sph, cyl to -4.00
Fitting Height	18mm
Width of Far Vision - 2mm above MRP	16.0mm
Minimum Width of Intermediate	7.0mm
Width of Near Zone - 14mm	11.5mm
Width of Near Zone - 20mm	15.2mm
Start of Progression above MRP	2.0mm
End of Progression below MRP	-14.0mm
Length of Progression	16.0mm



Base 5.25 Add 2.00